


AqT® Biotin Reduced Thiol-Labeled Trastuzumab (Anti-Her2 mAb)

Lyophilized Powder
100 µg or 500 µg of conjugate per tube,
lyophilized from a citrate buffer
containing a sugar-based stabilizer.
≥95% monomeric by SEC HPLC



Product Number:
CM86169

Product Description

Biotin is highly hydrophobic and only sparingly soluble in water. These properties make it difficult to achieve high biotin loading on a biopolymer. Most conventional biotinylation reagents require significant amounts of organic solvent, which can promote protein aggregation, precipitation, and loss of activity during labeling. CellMosaic's proprietary AqT® linker converts biotin into a highly water-soluble reagent, enabling high biotin loading with minimal or no organic solvent.

Trastuzumab is a humanized IgG1κ monoclonal antibody that binds to the extracellular domain of human epidermal growth factor receptor 2 (HER2/ErbB2). Trastuzumab is the antibody component of HER2-targeted antibody–drug conjugates such as ado-trastuzumab emtansine and trastuzumab deruxtecan. Biotinylated trastuzumab can be used as a research reagent for HER2 studies and diagnostic assay development. Potential biotin–avidin-based applications include ELISA, flow cytometry, immunohistochemistry (IHC), immunocytochemistry (ICC), immunoprecipitation (IP), pull-down assays, biosensor and surface-immobilization applications, cell isolation, and magnetic separation.

Unlike conventional biotin NHS ester chemistry, which randomly labels surface lysine residues, this AqT® biotinylated trastuzumab uses CellMosaic's proprietary super-hydrophilic AqT® linker to provide site-specific biotinylation at reduced Fc-region disulfide bonds. This site-specific labeling approach is designed to minimize disruption of antigen binding and increase streptavidin-mediated signals. AqT® biotinylated trastuzumab contains an average of 6–8 biotin molecules per antibody, providing enhanced signal amplification, improved assay sensitivity, and greater signal intensity across a range of applications.

The product is supplied in one vial containing either 100 µg (Cat. No. CM86169-100UG) or 500 µg (Cat. No. CM86169-500UG). For bulk orders, please contact us for a quote.

Product Information

- **Chemical Name:** AqT® Biotin Reduced Thiol-Labeled Trastuzumab (Anti-Her2 mAb)
- **Molecular Weight:** Approximately 152 kDa
- **Antibody Information:** Humanized monoclonal IgG1κ antibody targeting human HER2/ErbB2
- **Physical Appearance:** White to off-white lyophilized powder
- **Storage Temperature:** -20 °C or below
- **Purity:** ≥95% monomeric by SEC-HPLC (refer COA of each lot for actual value)
- **Average Biotin Loading:** 6–8 biotin molecules per trastuzumab molecule (refer COA of each lot for actual value)

Key Features

- Lyophilized powder that is ready for use after reconstitution with water; no additional buffer is required
- High biotin loading, with an average of 6–8 biotin molecules per trastuzumab molecule
- Protein content accurately determined by UV/HPLC analysis; biotin incorporation levels determined by biotin assay
- **Advanced AqT® linker technology:** Utilizes CellMosaic's proprietary, super-hydrophilic, water-soluble, charge-neutral AqT® linker to enhance aqueous solubility and minimize aggregation

Reconstitution

Reconstitute with the desired volume of sterile water. For reference, reconstituting one 100 µg tube with 100 µL yields a product concentration of 1 mg/mL. The resulting solution is citrate-buffered at approximately pH 6.4.

Applications for HER2

- HER2 detection and quantification
- HER2 binding and competition assays
- Selection and capture of HER2-positive cells
- Affinity purification and pull-down assays for HER2 and associated protein complexes
- HER2 receptor internalization studies
- Targeted imaging of HER2-positive cells and tumors
- HER2-targeted delivery research
- Assay development and quality control for HER2-expressing cell lines, biosensors, diagnostic assay-development platforms, and HER2-targeted formulations

Storage/Stability

- Store the lyophilized product at or below -20 °C
 - Shelf life before reconstitution: One year from the date of receipt when stored as directed
- After reconstitution, store at 2-8 °C and use within 24 hours for best performance

Figure 1. Representative size-exclusion HPLC (SEC-HPLC) profiles of trastuzumab before labeling (blue trace) and AqT® Biotin Trastuzumab after labeling (red trace; lot no. S702.S11.0701F; DOL: 6.9) in citrate buffer, demonstrating ≥99% monomeric purity with no detectable aggregates.





Figure 2. Representative hydrophobic interaction chromatography (HIC-HPLC) profiles of trastuzumab before labeling (blue trace) and AqT® Biotin Trastuzumab after labeling (red trace; lot no. S702.S11.0701F; DOL: 6.9) in citrate buffer, showing only a slight hydrophobic shift and a high degree of antibody homogeneity.



Important Notes

For Research Use Only. Not for Use in Human or Diagnostic Procedures.

The information contained in this document is subject to change without notice. CellMosaic assumes no responsibility for any errors or omissions in this document. Under no circumstances shall CellMosaic be liable, whether under contract, tort, warranty, statute, or any other legal theory, for any special, incidental, indirect, punitive, multiple, or consequential damages related to or arising from this document or the use of this product.

Intellectual Property

CellMosaic's products, including but not limited to AqT®, Nelon™, oxLink™ and sxLink™ (collectively, the "Products"), are protected by issued patents and pending patent applications in the United States and other jurisdictions. A current list of applicable patents is available upon request.

Restricted Use

The purchase of the Product conveys only a limited, non-exclusive, non-transferable, and non-sublicensable license to use the Product for the purchaser's internal research and development purposes. No right or license is granted, whether expressly or by implication, estoppel, or otherwise, under any patent, trademark, copyright, or other intellectual property right of CellMosaic, except as expressly stated herein.

Any commercial use of the Product is strictly prohibited without a separate written agreement with CellMosaic. Prohibited commercial uses include, without limitation: (a) the sale, lease, licensing, distribution, or other transfer of the Product or any materials derived from or produced using the Product; (b) the sale, lease, licensing, or other grant of rights to use the Product or any materials derived from or produced using the Product; and (c) the use of the Product to perform services for third parties for a fee or other consideration, including, without limitation, contract research, screening services, and diagnostic or therapeutic applications

For information regarding commercial licensing or collaboration, please contact us at info@cellmosaic.com.